

RoboCup@Home Practical course

Tutorials

Dr. Karinne Ramirez-Amaro

Dr. Emmanuel Dean

Dr. Pablo Lanillos


M.Sc. Roger Guadarrama

M.Sc. Constantin Uhde

Dr. Gordon Cheng

Technical University of Munich

Chair for Cognitive Systems



Definition of teams

- 1) **Hands-on tutorials:** introductory tutorials to get familiar with the equipment (robot & sensors) used for this course.
- 2) **Group definition tasks:** form groups according to the students' knowledge to address different challenges of the competition. The students designate a team leader.
- 3) **Development and test phase:** design and implementation of algorithms to solve the problems defined for your working team.
- 4) **Final phase:** test real scenarios on a mobile robot to evaluate the performance of the robots abilities.

Definition of teams

2) Group definition tasks: form groups according to the students' knowledge to address different challenges of the competition.

- Designate a team leader.

Responsibilities of the team leader:

- **Coordinate** the work of the group. Make sure that the work is correctly **distributed**.
- Be responsible for the **key** of the laboratory
- Direct **communication** with the supervisor(s)

Definition of topics

Category III: This category includes:

- following a human,
- indoor navigation in crowded environments,
- recognizing & grasping alike objects,
- find a calling person (waving or shouting), etc.
- deal with incomplete information

IMPORTANT: Define the strongest capabilities of your team:
Control, navigation, perception, and/or learning.

Tutorial 5: TIAGo – getting started

Tutorial 5: TIAGo – getting started

TIAGo overview



©Pictures taken from PAL Robotics documentation

TIAGo T.S.

Tutorial 5: TIAGo – getting started

Mobile base



Laser range-finder



Service panel:
HDMI, USB 3.0, On/Off computer



Rear sonars


Charging connector entry



©Pictures taken from PAL Robotics documentation

Tutorial 5: TIAGo – getting started

User panel




©Pictures taken from PAL Robotics documentation

Tutorial 5: TIAGo – getting started

Torso lifter



Stroke: 350 mm
Max speed: 50 mm/s




©Pictures taken from PAL Robotics documentation

TIAGo T.S.

Tutorial 5: TIAGo – getting started

Arm and end-effector



Force/torque sensor



Hey5 hand




PAL gripper



©Pictures taken from PAL Robotics documentation

Tutorial 5: TIAGo – getting started



Arm reach



©Pictures taken from PAL Robotics documentation

Tutorial 5: TIAGo – getting started

Head



Tutorial 5: TIAGo – getting started

Start the robot



Turn the electric key



Release the emergency button



Press On button during
1 second



©Pictures taken from PAL Robotics documentation

Tutorial 5: TIAGo – getting started

Charging TIAGo

- Use the charger supplied with TIAGo in order to charge the robot



Open the lid located
in the rear part



Insert charging
connector with metal lock
facing up, push until you
hear a 'click'



The duration of the charging
is about 4 h for 1 battery
completely discharged, and
about 8 h for 2 batteries



Insert charging
connector with metal
lock facing up, push
until you hear a
'click'



Tutorial 5: TIAGo – Safety measures

Safety measures




Tutorial 5: TIAGo – Safety measures

Proper robot shut down



Execute **Home** predefined motion through i.e. the Web commander



Lower the torso as much as possible using the joystick LT button



Hold the wrist of the robot ...

Keep holding the wrist until the motors are powered off. Then place carefully the wrist on top of the mobile base cover




... and press the On/Off button for 1 second

©Pictures taken from PAL Robotics documentation

TIAGo T.S.

Tutorial 5: TIAGo – Safety measures

How to get the arm out of Home safely



Get the arm out of collision:

<http://tiago-0c:8080>

5. Demos

Get out of collision

Raise the torso to its maximum height with the joystick



Execute the movement **unfold_arm** through i.e. the Web commander

<http://tiago-0c:8080>

6. Movements


unfold_arm



Tutorial 5: TIAGo – Safety measures

Measures to prevent the robot falling (I)

- The robot has been designed to be stable even when the arm is holding its maximum payload in its most extreme kinematic configuration
- Nevertheless, some measures need to be respected in order to avoid the robot falling



Do not apply external forces to the arm when it is extended and holding a weight



©Pictures taken from PAL Robotics documentation

Tutorial 5: TIAGo – Safety measures

Measures to prevent the robot falling (II)



Do not navigate at **high speeds** with the arm extended and holding a weight specially when the torso is extended



©Pictures taken from PAL Robotics documentation

Tutorial 5: TIAGo – Safety measures

Measures to prevent the robot falling (III)

- Do not navigate on floors with slope higher than 6°
- And in that case navigate with the torso at its lowest level and with the arm folded, i.e. in **Home** configuration



©Pictures taken from PAL Robotics documentation

Tutorial 5: TIAGo – Safety measures

Measures to achieve a safe navigation

- When carrying a weight:



It is highly recommended to navigate with the arm folded and the torso at low extension, like in the **Home** configuration




- Reduces the footprint of the robot and lowers the probability to collide with the environment with the arm
- Ensures the stability of the robot



©Pictures taken from PAL Robotics documentation

Tutorial 5: TIAGo – Safety measures

How to proceed when a collision occurs



©Pictures taken from PAL Robotics documentation

Tutorial 5: TIAGo – Software

Software architecture overview

- Operating system:

Ubuntu 14.04 LTS 64-bit



Xenomai real-time framework



real-time control
software

- Robotics middleware:

Orocos 2.8



ROS Indigo



PAL dubnium



ROS packages developed by PAL Robotics



Tutorial 5: TIAGo – Software

ROS packages

- **ROS** Indigo packages:

Installation path: **/opt/ros/indigo**

- **PAL** dubnium packages:

Installation path: **/opt/pal/dubnium**

- **User** packages:

Installation path: **/home/pal/deployed_ws**

Sub-folder	Description
bin	nodes (executables)
include	package header files
lib	package dynamic libraries and python files
share	packages cmake, launch and config files
etc	other files

All the ROS software in the robot is executed
with the user **pal**



Tutorial 5: TIAGo – Deploy new software in the robot

Connect to the robot ROS topics from the remote PC

On you local PC execute following commands:

1) `export ROS_MASTER_URI=http://tiago-24c:11311`


2) `export ROS_IP=ip_of_your_PC`

Tutorial 5: System integration on TIAGo

- **Finish all the tutorials of TIAGo**
- **Get really familiar with the Navigation package of TIAGo**

Tutorial 5: System integration

Exercise 1: Build a map of one of the floors from ICS. You can choose between the 2nd floor or the basement. The robot should be able to navigate within these rooms.



Tutorial 5: System integration

Exercise 2

- **Object Recognition**
- **Robot Manipulation**
- **Mapping and Navigation**

Tutorial 5: System integration

Exercise 2: The robot needs to move to the shelf grab one item and move it to the table.

Perception:


- 3D object localization
- Object detection
- Identify the shelf/table

Robot manipulation:

- Command the robot arm to a desired pose
- Grasp the object

Navigation:

- Move to the shelf and avoid collisions
- Move to the table





Tutorial 5: System integration

- **Object Recognition**


Perception

From plane segmentation to table computation




`pcl::ModelCoefficients::Ptr`

$$ax + by + cz + d = 0$$



Perception


Segmenting the objects on the bookshelf



- Transform to frame if needed (`pcl_ros::transformPointCloud`)
- Filter limits (`Passthrough`)
- Downsample (`VoxelGrid`)
- Remove the planes (`pcl::SACSegmentation`)
- Find the cylinders (`pcl::SACSegmentationFromNormals`)
- Remove outliers (`pcl::StatisticalOutlierRemoval`)
- Publish the list of pointclouds or the one needed

Perception

What else we can get from the Point Cloud data?



- Shape coefficients: cylinder [point_on_axis (x,y,z), axis_direction, (x,y,z) radius]
- 3D centroid
- Orientation
- Bounding box: min (x,y,z) and max (x,y,z) | centre (x0,y0,z0)
- Colored Texture: `VoxelGrid.setDownsampleAllData(true);`
- Descriptors
- Distance to the object

In order to grasp we need the **pose** of the object.

Perception

Tips for visualization in rviz:

- One object: send the PointCloud2
- Several objects: use markers (visualization_msgs/MarkerArray.msg)

```
#ifndef PUBLISH_MARKERS
visualization_msgs::Marker ObjectClustering::getCloudMarker(const pcl::PointCloud<PointType>::Ptr cloud, int id)
{
    //create the marker
    visualization_msgs::Marker marker;
    marker.header.frame_id = processing_frame_;
    marker.header.stamp = ros::Time();
    marker.action = visualization_msgs::Marker::ADD;
    marker.lifetime = ros::Duration(5);
    marker.ns = "segmentation";
    marker.id = id;
    marker.pose.orientation.w = 1;

    marker.type = visualization_msgs::Marker::POINTS;
    marker.scale.x = 0.002;
    marker.scale.y = 0.002;
    marker.scale.z = 1.0;

    marker.color.r = ((double)rand())/RAND_MAX;
    marker.color.g = ((double)rand())/RAND_MAX;
    marker.color.b = ((double)rand())/RAND_MAX;
    marker.color.a = 1.0;


    for(size_t i=0; i<cloud->size(); i++) {
        geometry_msgs::Point p;
        p.x = (*cloud)[i].x;
        p.y = (*cloud)[i].y;
        p.z = (*cloud)[i].z;
        marker.points.push_back(p);
    }

    return marker;
}
#endif
```

Tutorial 5: System integration

- **Object Recognition**
- **Robot Manipulation**

Movel! Motion planner



Grasping with MoveIt! planner

- 1) Define object position (From perception node)
- 2) Feed the object position and geometry restrictions to the move_group node.
- 3) Request a planning service.
- 4) Wait for response.
- 5) Execute the motion.
- 6) Check MoveIt! Tutorial from <http://wiki.ros.org/Robots/TIAGo/Tutorials>

Tutorial 5: System integration

- **Object Recognition**
- **Robot Manipulation**
- **Mapping and Navigation**

Navigation and Mapping


Hints:

- Generate a map for the testing environment.
- Use the obtained map to do navigation.
- Use MoveBaseGoal to move the robot with respect to the base_link.

If you want to navigate in the unknown environment (without previously building a map) you can check another SLAM packages :

http://wiki.ros.org/hector_slam or http://wiki.ros.org/slam_karto

Tutorial 5: System integration



Thank you

robocup.atHome.ics@tum.de

www.ics.ei.tum.de